

Dräger



This is Ventilation

Evita® V800 and Evita® V600

Improving outcomes in Intensive Care

We focus on reducing mortality rates in the ICU, supporting patient outcome and increasing staff satisfaction in the ICU via connected technologies and services that help achieve therapeutic goals faster and safer.



That's why we are
**Your Specialist
in Acute Care**

Your Specialist in Acute Care since 1889



1907 |
Pulmotor

The first life-saving ventilator



Today |
Evita V600/V800

Experience the next level of ventilator operation. The Evita® V600/V800 combines high performance ventilation with an aesthetic design enabling quick and efficient operation.

We share ...

Improving
clinical outcomes



Managing
cost of care



...
your goals
...

Ensuring
staff satisfaction



Enhancing patient
experience



Our mission – Improving outcomes in Intensive Care

As your specialist in acute care

we focus on reducing mortality rates in the ICU. Supporting patient outcome and increasing staff satisfaction in the ICU via connected technologies and services that help achieve therapeutic goals faster and safer is what drives us.


Only 41%

of patients who are artificially ventilated for at least 14 days, will survive the next year.¹

Avoid ICU-acquired weakness

Start the weaning process as early as possible to help reduce ventilation time, improve sedation management and optimise patient interaction.

Encourage spontaneous breathing, which helps train of respiratory muscles and start early mobilisation activities.



Prevention of VALI / ARDS

Protect the lung with personalised ventilation strategies. Support spontaneous breathing at any time to provide a seamless transition from controlled ventilation to patient contribution. Improve clinical results with decision supporting insights from multiple sources.

Avoiding cognitive impairment

Provide a comfortable and supportive environment that helps your patient feel calm and at ease. Turn the ICU into a healing environment and enable the patient to feel more comfortable in a family-friendly integrated surrounding.

01 _____

Operation & Handling

Operation & Handling

Clear information on the user interface are key to follow a defined strategy in ventilating a patient.

We support an intuitive operation of our devices which lowers education times and helps to avoid human errors.



Operation & Handling

Quick and safe to operate

Quick and safe to operate even in the most stressful situations due to intuitive menu access to both settings and your clinical data.



Fully recorded

All patient data, alarms and trends are fully recorded. Conveniently exported via USB interface.



Step-by-step guidance

Step-by-step guidance leads you through every procedure.




Multiple view configurations

Switch between multiple view configurations with the touch of a finger.

|

Individual protective ventilation therapy along The Respiration Pathway



As non-invasive as possible, as invasive as necessary: The variability and diversity of our treatment tools for use along the respiration pathway enable you to administer protective mechanical ventilation therapy in your ICU.

02

Safeguarding High-Flow O₂-Therapy

Evita Ventilators – Safeguarding O₂-Therapy

Safeguarding

Our O₂-Therapy function enables a Safeguarding High Flow Therapy as you have the possibility of a **Pmax setting**. A possible exceeding of the set maximum allowed airway pressure is detected early, and **the flow is reduced accordingly**.

– For all patients


O₂-Therapy (high flow oxygen therapy) is suitable for use as respiratory support for adults, paediatric patients, and neonates who can inhale and exhale spontaneously.

– Seamless transition

By using Dräger ventilators you can conveniently switch your patient between traditional invasive ventilation, non-invasive ventilatory support or HFNC to meet changing needs while still using the same device.



Life-Saving Treatment...



...with High-Flow
Oxygen Therapy
for All Patient
Categories.




The Respiration Pathway

Hi-Flow Star Cannula

The Respiration Pathway

Support whenever you need it



Prevent

Avoiding intubation as long as possible helps to decrease the possibility of patients developing VAP.

Stabilise

Stabilising your patient is crucial to protecting their lungs from potential conditions such as VALI, which could result in ALI or ARDS.

Wean

Weaning your patient efficiently reduces incidences of ALI or ARDS.

Recover

Ensuring successful weaning and avoiding reintubation improves your patient's recovery.

Respiration Pathway

As non-invasive as possible, as invasive as necessary:

Discover the variability and diversity of our treatment tools for use along the respiration pathway. They enable you to administer protective mechanical ventilation therapy in your ICU.

Physiology and Physics of High Flow Oxygen Therapy

Physiology, HFOT provides a consistent and high concentration of warmed, humidified oxygen that meets or exceeds patients' inspiratory flow demands, reducing work of breathing, decreasing dead space, and improving oxygenation by increasing alveolar ventilation.

From a **physics** perspective, HFOT utilizes principles such as laminar flow to deliver a high and stable flow rate of gas through large-bore nasal cannulas, creating a slight positive airway pressure that helps keep alveoli open, reducing atelectasis and improving gas exchange.

The high flow also flushes carbon dioxide from the upper airway, decreasing rebreathing and dead space, which enhances overall ventilation efficiency.




Literature – FLORALI & Frat et al., 2015

The FLORALI trial ([Fluid and Oxygen in Patients with Acute Lung Injury](#)) primarily investigated the effectiveness of high-flow nasal cannula (HFNC) oxygen therapy compared to standard oxygen therapy and non-invasive ventilation in patients with acute hypoxemic respiratory failure.

It found that HFNC was associated with [lower intubation rates](#) and improved patient outcomes, especially in [cases of moderate hypoxemia](#).

The study by Frat et al. published in the New England Journal of Medicine 2015. The trial demonstrated high-flow oxygen therapy's effectiveness in reducing intubation rates among patients with hypoxemic respiratory failure, specifically in the context of acute respiratory failure due to [pneumonia](#), [cardiogenic pulmonary edema](#), or other causes.



Literature – ROX-Index & Christophe Girault et al., 2024

Recent reviews from 2023–24 emphasize that flow rates in respiratory support devices, such as HFNC, significantly influence treatment effectiveness, with higher flows improving oxygenation and reducing work of breathing.

Christophe Girault et al. published their study on the performance of the **ROX index** in predicting high-flow nasal oxygen outcomes in COVID-19-related hypoxemic acute respiratory failure in the journal Respiratory Medicine.

➤ [Ann Intensive Care](#). 2024 Jan 18;14(1):13. doi: 10.1186/s13613-023-01226-6.

ROX index performance to predict high-flow nasal oxygen outcome in Covid-19 related hypoxemic acute respiratory failure

Christophe Girault ^{1 2}, Michael Bubenheim ³, Déborah Boyer ⁴, Pierre-Louis Declercq ⁵, Guillaume Schnell ⁶, Philippe Gouin ⁷, Jean-Baptiste Michot ⁸, Dorothée Carpentier ⁴, Steven Grangé ⁴, Gaëtan Béduneau ⁹, Fabienne Tamion ¹⁰

Affiliations + expand

PMID: 38236356 PMID: [PMC10796865](#) DOI: [10.1186/s13613-023-01226-6](#)

Abstract

Background: Given the pathophysiology of hypoxemia in patients with Covid-19 acute respiratory failure (ARF), it seemed necessary to evaluate whether ROX index (ratio $\text{SpO}_2/\text{FiO}_2$ to respiratory rate) could accurately predict intubation or death in these patients initially treated by high-flow nasal oxygenation (HFNO). We aimed, therefore, to assess the accuracy of ROX index to discriminate

Downloaded by 90.196.33.117 from www.liebertpub.com at 09/23/25. For personal use only.



Simon Oczkowski^{1,2,26}, Begüm Ergen^{3,26}, Lieuwe Bos^{4,5}, Michelle Chatwin⁶, Miguel Ferrer⁷,
Cesare Gregoret^{8,9}, Leo Heunks¹⁰, Jean-Pierre Frat^{11,12}, Federico Longhini¹³, Stefano Nava^{14,15},
Paolo Navale^{16,17}, Aylin Özşancak Uğurlu¹⁸, Lara Pisanì^{14,15}, Teresa Renda¹⁹, Arnaud W. Thille²⁰,
Wolfram Windisch²¹, Thomý Tonio²², Jeanette Boyd²³, Giovanni Sotgiu²⁴

Frederic Vargas MD PhD, Mélanie Saint-Leger MD, Alexandre Boyer MD PhD,
Nam H Bui MD, and Gilles Hilbert MD PhD

INTRODUCTION: High-flow nasal cannula (HFNC) can deliver heated and humidified gas (up to 100% oxygen) at a maximum flow of 60 L/min via nasal prongs or cannula. The aim of this study was to assess the short-term physiologic effects of HFNC. Inspiratory muscle effort, gas exchange, dyspnea score, and comfort were evaluated. **METHODS:** Twelve subjects admitted to the ICU for acute hypoxic respiratory failure were prospectively included. Four study sessions were performed. The first session consisted of oxygen therapy given through a high- F_{IO_2} -non-rebreathing face mask. Recordings were then obtained during periods of HFNC and CPAP at 5 cm H_2O in random order, and final measurements were performed during oxygen therapy delivered via a face mask. Each of these 4 periods lasted ~20 min. **RESULTS:** Esophageal pressure signals, breathing pattern, gas exchange, comfort, and dyspnea were measured. Compared with the first session, HFNC reduced inspiratory effort (pressure-time product of $156.0 [119.2-194.4]$ cm H_2O $\times$ s/min vs $204.2 [149.6-324.7]$ cm H_2O $\times$ s/min, $P < .01$) and breathing frequency ($P < .01$). No significant differences were observed between HFNC and CPAP for inspiratory effort and breathing frequency. Compared with the first session, P_{aO_2}/F_{IO_2} increased significantly with HFNC ($167 [157-184]$ mm Hg vs $156 [110-171]$ mm Hg, $P < .01$). CPAP produced significantly greater P_{aO_2}/F_{IO_2} improvement than did HFNC. Dyspnea improved with HFNC and CPAP, but this improvement was not significant. Subject comfort was not different across the 4 sessions. **CONCLUSIONS:** Compared with conventional oxygen therapy, HFNC improved inspiratory effort and oxygenation. In subjects with acute hypoxic respiratory failure, HFNC is an alternative to conventional oxygen therapy. (ClinicalTrials.gov registration NCT01056952) **Key words:** high-flow nasal cannula; continuous positive airway pressure; oxygen therapy; acute hypoxic respiratory failure; inspiratory effort. [Respir Care 2015;60(10):1369-1376. © 2015 Daedalus Enterprises]

of Critical Care, McMaster University, Hamilton, ON, Canada. ²Dept of Health Research Methods, Hamilton, ON, Canada. ³Dept of Pulmonary and Critical Care, Dokuz Eylul University, Dept of Intensive Care and Laboratory of Experimental Intensive Care and Anesthesiology, Academic Medical Center, Amsterdam, The Netherlands. ⁴Respiratory Medicine, Amsterdam University Medical Center, The Netherlands. ⁵Academic and Clinical Department of Sleep and Breathing Research Unit, Royal Brompton and Harefield NHS Foundation Trust, London, UK. ⁶Dept of Pulmonary and Critical Care, IDIBAPS, University of Barcelona and CIBERes, Barcelona, Spain. ⁷Dept of Surgical Sciences, University of Palermo, Italy. ⁸G. Giglio Institute, Cefalù, Italy. ⁹Dept of Intensive Care Medicine, Amsterdam, The Netherlands. ¹⁰Centre Hospitalier Universitaire de Poitiers, Médecine Intensive Réanimatoire, Centre d'Investigation Clinique 1402 ALIVE, Université de Poitiers, Poitiers, France. ¹¹Anesthesiology and Surgical Sciences, Magna Graecia University, Catanzaro, Italy. ¹²Dept of Clinical Medicine, Alma Mater Studiorum University of Bologna, Bologna, Italy. ¹³IRCCS Azienda Ospedaliera di Ricovero e Cura, Respiratory and Critical Care Unit, Bologna, Italy. ¹⁴Department of Medicine – DIMED, University of Padua, Padua, Italy. ¹⁵Dept of Pulmonary Medicine, Padua University Hospital, Padua, Italy. ¹⁶Dept of Pulmonary Medicine, Padua University Hospital, Padua, Italy. ¹⁷Cardiothoracic and Vascular Department, Respiratory and Critical Care Unit, Careggi University Hospital, Florence, Italy. ¹⁸Department of Medicine, Universidade do Porto, Porto, Portugal. ¹⁹Dept of Pneumology, Cologne Merheim Hospital, Cologne, Germany. ²⁰Department of Medicine, Cologne Merheim Hospital, Cologne, Germany. ²¹Department of Medicine, Witten/Herdecke University, Faculty of Health/School of Medicine, Köln, Germany. ²²Department of Medicine, University of Bern, Bern, Switzerland. ²³European Lung Foundation (ELF), Sheffield, UK. ²⁴Statistics Unit, Dept of Medical, Surgical, Experimental Sciences, University of Sassari, Sassari, Italy. ²⁵Intensive Care Unit, Cardio-Thoraco-Neuro-Vascular Dept, Usl Toscana Sudest, S. Donato, Italy. ²⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ²⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ²⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ²⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ³⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁴⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁵⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁶⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁷⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁸⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹¹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹²Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹³Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁴Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁵Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁶Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁷Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁸Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ⁹⁹Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy. ¹⁰⁰Intensive Care Unit, Usl Toscana Sudest, S. Donato, Italy.

Corresponding author: Raffaele Scala (raffaele_scala@hpi.de)

Shareable abstract (@ERSpublications)

This guideline provides evidence-based recommendations for the use of high-flow nasal cannula alongside other noninvasive forms of respiratory support in adults with acute respiratory failure.

Cite this article as: Oczkowski S, Ergan B, Bos L, et al. ERS clinical practice guidelines: high-flow nasal cannula in acute respiratory failure. *Eur Respir J* 2022; 59: 2101574 [DOI: 10.1183/13993003.000000002101574].

British Journal of Anaesthesia 103 (6): 886–90 (2009)
doi:10.1093/bja/aep280 Advance Access publication October 20, 2009

BJA

Nasal high-flow therapy delivers low level positive airway pressure

R. Parke^{1*}, S. McGuinness¹ and M. Eccleston^{2,†}

¹Cardiothoracic and Vascular Intensive Care Unit, Level 4, Building 32, Auckland City Hospital, Private Bag 92024, Auckland 1010, New Zealand. ²Fisher and Paykel Healthcare, 1 Maurice Paykel Place, East Tamaki, Auckland, New Zealand

*Corresponding author. E-mail: rparke@adhb.govt.nz

Background. The aim of this prospective study was to determine whether a level of positive airway pressure was generated in participants receiving nasal high flow (NHF) delivered by the Optiflow™ system (Fisher and Paykel Healthcare Ltd, Auckland, New Zealand) in a cardiothoracic and vascular intensive care unit (ICU).

Methods. Nasopharyngeal airway pressure was measured in 15 postoperative cardiac surgery adult patients who received both NHF and standard facemask therapy at a flow rate of 35 litre min⁻¹. Measurements were repeated in the open mouth and closed mouth positions. Mean airway pressure was determined by averaging the pressures at the peak of inspiration of each breath within a 1 min period, allowing the entire pressure profile of each breath to be included within the calculation.

Results. Low level positive pressure was demonstrated with NHF at 35 litre min⁻¹ with mouth closed when compared with a facemask. NHF generated a mean nasopharyngeal airway pressure of mean (SD) 2.7 (1.04) cm H₂O with the mouth closed. Airway pressure was significantly higher when breathing with mouth closed compared with mouth open ($P<0.0001$).

Conclusions. This study demonstrated that a low level of positive pressure was generated with NHF at 35 litre min^{-1} of gas flow. This is consistent with results obtained in healthy volunteers.

Australian Clinical Trials Registry www.actr.org.au ACTRN012606000139572

Br J Anaesth 2009; 103: 886–90

Keywords: airway pressure; nasal high-flow therapy; Optiflow™; oxygen


Accepted for publication: August 18, 2009

Patients with respiratory failure are typically treated with three main respiratory support strategies, depending on the severity of their illness. These are traditional oxygen therapy, non-invasive ventilation, and invasive mechanical ventilation.¹ A new respiratory support therapy has recently been introduced into the adult arena. NHF allows the delivery of up to 60 litre min⁻¹ of heated and humidified gas via a wide bore nasal cannula. However, the

been reports that NHF may be beneficial in the treatment of obstructive sleep apnoea, attributable to a flow-related pressure effect.⁷ There has been a similar pressure effect reported in healthy adult volunteers⁸⁻⁹ where a positive relationship between flow and airway pressure has also been described. However, to date there have been no

High-Flow Oxygen Therapy – Mechanisms of action

- Flushing of anatomical dead space
- Reduced work of breathing
- Enhanced humidification:
 - Improved mechanisms
 - Reduction in the metabolic cost
 - Better mucociliary clearance
- Better control of the patient's FiO_2
- Positive airway pressure improves alveolar recruitment



An Experimental and Numerical Investigation of CO₂ Distribution in the Upper Airways During Nasal High Flow Therapy

S. C. VAN HOVE,¹ J. STOREY,¹ C. ADAMS,² K. DEY,² P. H. GEOGHEGAN,² N. KABALIUK,² S. D. OLDFIELD,³
C. J. T. SPENCE,³ M. C. JERMY,² V. SURESH,^{1,4} and J. E. CATER¹

¹Department of Engineering Science, University of Auckland, Auckland 1142, New Zealand; ²Department of Mechanical Engineering, University of Canterbury, Christchurch 8041, New Zealand; ³Fisher and Paykel Healthcare Limited, 15 Maurice Paykel Place, Auckland 2013, New Zealand; and ⁴Auckland Bioengineering Institute, University of Auckland, Auckland 1142, New Zealand

(Received 14 December 2015; accepted 31 March 2016; published online 8 April 2016)

Associate Editor Kerry Hourigan oversaw the review of this article.

Combines computational modelling and experimental measurements to analyse how nasal high-flow (NHF) therapy influences CO₂ clearance in the upper respiratory tract.

An Experimental and Numerical Investigation of CO₂ Distribution in the Upper Airways During Nasal High Flow Therapy

S. C. VAN HOVE,¹ J. STOREY,¹ C. ADAMS,² K. DEY,² P. H. GEOGHEGAN,² N. KABALIUK,² S. D. OLDFIELD,³
C. J. T. SPENCE,³ M. C. JERMY,² V. SURESH,^{1,4} and J. E. CATER¹

¹Department of Engineering Science, University of Auckland, Auckland 1142, New Zealand; ²Department of Mechanical Engineering, University of Canterbury, Christchurch 8041, New Zealand; ³Fisher and Paykel Healthcare Limited, 15 Maurice Paykel Place, Auckland 2013, New Zealand; and ⁴Auckland Bioengineering Institute, University of Auckland, Auckland 1142, New Zealand

(Received 14 December 2015; accepted 31 March 2016; published online 8 April 2016)

Associate Editor Kerry Hourigan oversaw the review of this article.

The research demonstrates that NHF improves CO₂ removal by flushing out exhaled gases from the upper airways, thereby reducing dead space and enhancing ventilation efficiency. This mechanistic insight supports the physiological benefits of NHF in promoting better gas exchange, highlighting its role in reducing hypercapnia and improving patient outcomes during respiratory support, as well as informing optimal flow settings for effective therapy.

Velocity Distribution Over Sagittal Plane


Top – 0 L/M

Middle – 30 L/M

Bottom – 60 L/M

Left - Midway through expiration

Right – Instant in time prior to the transition between expiration and inspiration



CO₂ Volume Fraction Displayed Over Sagittal Plane







Top – 0 L/M

Middle – 30 L/M

Bottom – 60 L/M

Left - Midway through expiration


Right – Instant in time prior to the transition between
expiration and inspiration



Conclusion


“Numerical simulation of the flow and CO₂ concentration within the upper airway of a subject undergoing nasal high flow therapy have been performed. The numerical simulations were validated using an experimental model and the results are consistent with previously reported airway measurements.

A nasopharyngeal washout phenomenon was observed and the amount of CO₂ cleared from the nasal cavity was found to increase in a flow dependant manner over the range of flows simulated.” (Van Hove et al., 2016.)



Hi-Flow Star Nasal
Cannula: Straight
Prong Design

Hi-Flow Star Sizes



Hi-Flow Star Nasal Cannula: Correct Size

Use the sizing gauge to find the appropriate size.



Make sure that the prongs only **cover max. 50%** of the nostrils.



Hi-Flow Star Easy Application



Many thanks

Stephen Abram | Clinical Application Specialist